RIKTA

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DOUCHE 2

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DEAR CUSTOMER!

We congratulate you on purchasing the
EMITTER DOUCHE 2
and we express confidence that it will help
to maintain your health for many years!

The emitter is manufactured at the enterprise of JSC „MILTA – PCP GIT“, whose
quality management system complies with the requirements of international
standards.

The emitter is certified in the Russian Federation jointly with RIKTA®
device. Certificate of conformity № РОСС RU.ИМ02.Н17823.

Validity — until 26.08.2019 year.

**CONTENTS**

RU	3
EN	
1. GENERAL INFORMATION	12
2. INTENDED USE.....	13
3. DELIVERY SET	13
4. SAFETY PRECAUTIONS	14
5. DEVICE.....	15
6. TECHNICAL SPECIFICATION	16
7. OPERATING INSTRUCTION	17
8. STORAGE AND TRANSPORTATION	17
9. MANUFACTURER’S WARRANTY	18
FR	19
DE	27
WARRANTY CARD	35
WARRANTY REPAIRS	35

1. GENERAL INFORMATION

- The present operating manual (further – the Manual) is the document that certifies basic parameters and technical specifications of the emitters of models DOUCHE 2-02, DOUCHE 2-03, DOUCHE 2-04, DOUCHE 2-05 (further emitter DOUCHE 2), guaranteed by the Manufacturer, and is intended for acquaintance with the operating rules of the Emitter.

The Operating manual should be kept during all the Emitter service life.

- It is necessary to read the Manual carefully before using the emitter. Particular attention should be paid to compliance with the safety precautions listed on page 14.
- The methods of carrying out treatment procedures are stated in the Methodical manual on application of the emitter DOUCHE 2 for RIKTA® devices (further – the Methodical manual) or in individual doctor's prescription.
- The manufacturer reserves the right to make adjustments and amendments to this Manual and technical documentation at any time without notice.
- The manufacturer reserves the right to make changes in the configuration and design of the emitter, which will not affect its operation and basic technical specifications.

2. INTENDED USE

- The Emitter is intended for therapeutic noninvasive (without breaking the skin integument) treatment of scalp diseases.
- The Emitter should be used jointly with the magnetic-infrared laser therapeutic devices of RIKTA® series. Depending on the modification, the Emitter is intended for connection only to a specific model of the magnetic-infrared laser therapeutic device of RIKTA® series (the modification is specified in the nameplate on the emitter case):

DOUCHE 2-02 – with RIKTA®-02/1
(«QUANTERRA») device;

DOUCHE 2-03 – with RIKTA®-03/2 device;

DOUCHE 2-04 – with RIKTA®-04/4 device;

DOUCHE 2-05 – with RIKTA®-05 device.

3. DELIVERY SET

Delivery set of the Emitter DOUCHE 2

Description	Quantity, pcs
Emitter Douche 2	1
Operating manual	1
Methodical manual	1
Package	1

EN

4. SAFETY PRECAUTIONS

Using the Emitter is allowed only after acquaintance with this Manual and the Methodical manual.

- The Emitter is intended to work jointly with devices of RIKTA® series (see page 13). It is not acceptable to connect the Emitter to a device for which it is not intended!
- Don't open up the Emitter and don't attempt to disassemble or modify its internal components.
- During operation, it is recommended to avoid the impact of the direct or specularly reflected radiation into the eyes.

Don't stare into the beam and don't conduct the direct observations using optical instruments.

- Before and after the treatment procedure, wipe the Emitter with gauze or cotton, moistened in a disinfectant for medical devices based on 3% hydrogen peroxide solution.

Avoid any contact of the device with alcohol-containing means! Disinfection of the Emitter should be carried out only after switching off and disconnecting the device from the power supply.


The full list of safety precautions is specified in the Operating manual for your RIKTA® device. Indications and contraindications to laser therapy treatment are stated in the Methodical manual.

EN

The Emitter disposal shall meet the conventional standards and local regulations for the environment protection.

5. DEVICE

EN

6. TECHNICAL SPECIFICATION

Radiation wavelength, <i>mkm</i> :	
– laser infrared radiation.....	0,800...0,910
– broadband infrared radiation.....	0,860...0,960
– visible red radiation.....	0,600...0,700
Laser radiation pulse power, <i>W</i>	10^{+4}_{-2}
Infrared radiation average power, <i>mW</i>	60 ± 30
Red radiation average power, <i>mW</i>	7^{+3}_{-5}
Magnetic induction, <i>mT</i>	6 ± 2
Laser safety by GOST IEC 60825-1.....	class1M
Electrical safety by GOST P50267.0-92 / IEC 60601-1.....	type VP
Service life of the Emitter, years.....	5

7. OPERATING INSTRUCTION

Before connecting, make sure that the Emitter corresponds to your model of RIKTA® device (see page 13).

- Connect the cable [4] of the Emitter to the appropriate socket of the electronic control unit of your RIKTA® device.
- Turn on the device.

The operating instructions are described in the Operating manual of your RIKTA® device.

8. STORAGE AND TRANSPORTATION

- Each Emitter is packed in the Manufacturer's package.
- The Emitter should be stored in the Manufacturer's package in the heated room at temperature in the range of +10...+35°C and relative humidity of no more than 80% (at +25°C). The storage room should be free of dust, mercury vapor, acids and alkalis provoking corrosion.
- The Emitter should be transported by all types of the closed transport in Manufacturer's package. It is necessary to ensure a stable position of the package during transportation.

EN

9. MANUFACTURER'S WARRANTY

- The Manufacturer guarantees that the Emitter parameters correspond to the values specified in the Manual, within 24 months from the date of sale of the Emitter to the consumer (if the consumer complies with the conditions of use, transportation and storage). Within the aforesaid period, the manufacturer provides free repair of the device.
- The manufacturer is not responsible for and does not make free repairs of malfunctions caused by the fault of the consumer (for example, due to careless handling, transportation, non-compliance with the rules of use, as well as in case of breaking of the seal).
- This warranty does not apply to any of the following:
 - absence of the present *Manual*;
 - absence of the *Warranty card*;
 - mechanical damage including the damage during transportation;
 - breaking of the seal.
- The warranty also becomes void if during the warranty period repair of the faulty Emitter was carried out by unauthorized persons and organizations.
- After the warranty period, the manufacturer shall repair the Emitter at the expense of the consumer.
- In the absence of the sales date and the seller's stamp in the *Warranty card*, the warranty period shall be calculated from the date of manufacturing of the Emitter.
- All issues related to repairs and maintenance should be directed to the manufacturer or to the technical centers listed on the website www.rikta.ru or in a separate warranty booklet (if available).

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RIKTAMED LTD